



Laerdal

helping save lives



Airway management and CPR

Laerdal Advanced Airway Solution

The Laerdal Advanced Airway Solution is a fully equipped training platform for contaminated, basic, difficult airway management and CPR practice. Designed for engaging individual and team training, it delivers realism, real-time feedback, and fast setup and clean up. With adjustable airway complexity, it supports progressive learning by adapting to each learner's skill level and training goals.

Laerdal Advanced Airway Solution


Product Features

- ✓ Contaminated airway management with fluids e.g. simulated vomit
- ✓ Adjustable airway complexity features to match training outcomes
- ✓ Ready to use with all necessary training components included
- ✓ Easy and quick setup and clean-up under 5 minutes
- ✓ Realistic airway
- ✓ Suitable for individual and team training

Clinical skills

The Laerdal Advanced Airway Solution is a complete training solution designed to help learners practice everything from basic airway management and ventilation to advanced suction techniques, airway decontamination, and team based resuscitation.

Clinical skills include: contaminated airway, laryngospasm, pharyngeal swelling, tongue edema, trismus, head tilt - chin lift, lung aspiration, dynamic L/R lung resistance, chest compressions, bag mask ventilation, OPA, NPA



What's included

- Little Anne Advanced Airway Manikin (2)
- Contaminated Airway System, inc. not-for-clinical-use LSU (1)
- Thomas Select Tube Holder Adult (2)
- StiffNeck Select Extrication Collar (2)
- RQI Training BVM Adult (2)
- Reusable Laerdal Training Canister (2)
- Hand Pump System (2)
- Concentrated Simulated Vomit 75g
- Blood, CONC-RED 3oz
- Cleaning Adapter
- Cleaning cloth (2)
- Battery 1.5 V Alkaline (4)
- Airway Lubricant Glycerol (180ml)
- Silicone lubricant
- Inflation bulb assembly (2)
- Aspiration canisters (4)
- Tube clamp (2)
- Restraining strap (2)
- Counterweight container (4)
- Waterproof training mat large
- Organization pouch (4)
- Transport suitcase and strap

Ordering Information

13503150 Laerdal Advanced Airway Solution

The consumables kit includes:

252800 Concentrated simulated vomit x 2
252090 Airway Lubricant Glycerol
300-00750 Blood, CONC-RED
410-00050 Silicone lubricant

Advanced Protection Plan

135-000PP Advanced Protection Plan

