



Laerdal

helping save lives



Resusci® Anne Q CPR

More Efficient CPR Training

Resusci Anne Q CPR combines CPR and defibrillator training, airway management, and team-based resuscitation training in one integrated solution. With objective performance feedback through the Q CPR App, you can deliver realistic, measurable, guideline-compliant training that prepares learners for real cardiac arrest situations.

Resusci® Anne QCPR

In addition to providing high-quality CPR training, the Resusci Anne QCPR manikin allows for those skills and others to be integrated into team training. The available add-on options and accessories make Resusci Anne QCPR more than just a high-quality CPR training manikin. The availability of diverse skin tones, CPR bralette, QCPR app and integrated LINK technology at no added cost enable guidelines compliant and representative training.

Key Training Features

CPR

- Realistic compressions with adjustable resistance
- Ventilation training with volume feedback and chest rise

Defibrillation

- Automatic pad placement with ShockLink™ and AED Trainer
- Add your clinical AED/defibrillator using ShockLink™
- Realistic training with AED Trainer or Monitor by Laerdal

General

- Bluetooth Low Energy
- Rechargeable battery

QCPR App Feedback

The QCPR App transforms CPR training with real-time feedback, automatic result tracking, and easy performance comparisons at no added cost. Available on both iOS and Android devices, it supports a convenient bring-your-own-device approach.

- Compression depth
- Compression rate
- Hand placement
- Longest pause
- Ventilation volume
- Overall CPR score
- Progress over time



Clinical Features

	Standard Head	Airway Head
Airway occlusion	✓	✓
Jaw thrust Chin lift	✓	✓
Positive pressure ventilation	✓	✓
Head rotation for foreign airway obstruction	✓	✓
Mouth-to-mouth ventilation	✓	✗
Pocket mask ventilation	✓	✓
Bag-Valve-Mask ventilation	✓	✓
Manual Carotid Pulses	✓	✓
Oropharyngeal airway (OPA) placement	✗	✓
Nasopharyngeal airway (NPA) placement	✗	✓
Supraglottic airway (SGA) insertion	✗	✓
Endotracheal (ET) intubation	✗	✓
Anatomically accurate airway to the glottis	✗	✓
Sellick manoeuvre (cricoid pressure)	✗	✓
Stomach inflation simulation	✗	✓
Adjustable pupils	✗	✓

Ordering Information

177-06350 Resusci Anne QCPR Torso
Light Skin Tone

177-05350 Resusci Anne QCPR Full Body
Light Skin Tone

177-08350 Resusci Anne QCPR w/ Airway
Head – Torso - Light Skin Tone

177-07350 Resusci Anne QCPR w/ Airway
Head – Full Body - Light Skin Tone

177-06250 Resusci Anne QCPR Torso
Dark Skin Tone

177-05250 Resusci Anne QCPR Full Body
Dark Skin Tone

177-08250 Resusci Anne QCPR w/Airway
Head - Torso - Dark Skin Tone

177-07250 Resusci Anne QCPR w/ Airway
Head - Full Body - Dark Skin Tone

**Please contact the local sales team for medium skin tone orders.*

Optional Limb Accessories in Light Skin Tone

IV Arm, Trauma Limbs, Rescue Limbs, Bleeding Control Legs

Recommended Solutions

Monitor by Laerdal, AED Trainer and ShockLink

Feedback Devices

206-30001 SimPad® PLUS with SkillReporter

206-50050 SkillReporter Software License

170-30050 SkillGuide™

Warranty

One-year limited warranty. Refer to Laerdal Global Warranty for terms and conditions.

Laerdal.com



Laerdal
helping save lives